


Track and treat with confidence.



Flexor[®] Ansel

GUIDING SHEATH

Visibility

Radiopaque band, seamlessly incorporated within the sheath material, identifies precise location of the distal tip for positioning accuracy.

Trackability

Select configurations are available with hydrophilic coating to reduce friction and enhance trackability.

Kink-resistance

Patented coil-reinforced construction provides flexibility and resistance to kinking.




Flexor[®] Ansel

GUIDING SHEATH

Flexor[®] High Flex Ansel


Flexible enough for tortuous anatomy or crossover procedures, yet provides extra support when needed.


- Three-step construction enables flexibility for tight anatomy
- Stiffer shaft provides pushability
- Softer dilator improves trackability
- Check-Flo valve (4-8 Fr) or large Check-Flo valve (10-12 Fr)



Configurations

ANL0 


ANL1 

ANL2 

ANL3 

Flexor[®] Ansel

- Four different tip configurations to meet various anatomical needs
- Two-step construction provides flexibility
- Soft distal tip is atraumatic to the vessel
- Available with Tuohy-Borst Adapter or Check-Flo valve



Flexor[®] Ansel Guiding Sheath with Tuohy-Borst Side-Arm Adapter

Includes .018" and .038" compatible dilators

Tuohy-Borst

Order Number	Reference Part Number	Introducer Fr	Sheath ID inch/mm	Accepts Wire Guide Diameter inch	Sheath Length cm	Tip Configuration
G29995	KSAW-5.0-18/38-55-RB-ANL1-HC	5.0	.074/1.9	.018/.038	55	ANL1
G30120	KSAW-5.0-18/38-55-RB-ANL2-HC	5.0	.074/1.9	.018/.038	55	ANL2
G29996	KSAW-6.0-18/38-55-RB-ANL1-HC	6.0	.087/2.2	.038	55	ANL1
G29997	KSAW-6.0-18/38-55-RB-ANL2-HC	6.0	.087/2.2	.038	55	ANL2
G29999	KSAW-7.0-18/38-55-RB-ANL1-HC	7.0	.100/2.5	.038	55	ANL1
G30000	KSAW-7.0-18/38-55-RB-ANL2-HC	7.0	.100/2.5	.038	55	ANL2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Flexor® Ansel Guiding Sheath with Check-Flo® Hemostasis Valve and High-Flex Dilator

Includes .035" compatible dilator

Check-Flo Valve

Order Number	Reference Part Number	Introducer Fr	Sheath ID inch/mm	Accepts Wire Guide Diameter inch	Sheath Length cm	Tip Configuration
G56229	KCFW-4.0-35-55-RB-HFANL0-HC	4.0	.060/1.5	.035	55	ANL0
G56234	KCFW-4.0-35-55-RB-HFANL1-HC	4.0	.060/1.5	.035	55	ANL1
G56219	KCFW-4.0-35-70-RB-HFANL0-HC	4.0	.060/1.5	.035	70	ANL0
G56230	KCFW-5.0-35-55-RB-HFANL0-HC	5.0	.074/1.9	.035	55	ANL0
G56235	KCFW-5.0-35-55-RB-HFANL1-HC	5.0	.074/1.9	.035	55	ANL1
G56220	KCFW-5.0-35-70-RB-HFANL0-HC	5.0	.074/1.9	.035	70	ANL0
G56225	KCFW-5.0-35-70-RB-HFANL1-HC	5.0	.074/1.9	.035	70	ANL1
G56231	KCFW-6.0-35-55-RB-HFANL0-HC	6.0	.087/2.2	.035	55	ANL0
G56236	KCFW-6.0-35-55-RB-HFANL1-HC	6.0	.087/2.2	.035	55	ANL1
G56221	KCFW-6.0-35-70-RB-HFANL0-HC	6.0	.087/2.2	.035	70	ANL0
G56226	KCFW-6.0-35-70-RB-HFANL1-HC	6.0	.087/2.2	.035	70	ANL1
G56232	KCFW-7.0-35-55-RB-HFANL0-HC	7.0	.100/2.5	.035	55	ANL0
G49043	KCFW-7.0-35-55-RB-HFANL1-HC	7.0	.100/2.5	.035	55	ANL1
G56222	KCFW-7.0-35-70-RB-HFANL0-HC	7.0	.100/2.5	.035	70	ANL0
G56227	KCFW-7.0-35-70-RB-HFANL1-HC	7.0	.100/2.5	.035	70	ANL1
G56233	KCFW-8.0-35-55-RB-HFANL0-HC	8.0	.113/2.9	.035	55	ANL0
G49045	KCFW-8.0-35-55-RB-HFANL1-HC	8.0	.113/2.9	.035	55	ANL1
G56223	KCFW-8.0-35-70-RB-HFANL0-HC	8.0	.113/2.9	.035	70	ANL0
G56228	KCFW-8.0-35-70-RB-HFANL1-HC	8.0	.113/2.9	.035	70	ANL1
G56278	KCFW-10.0-35-45-RB-HFANL0-HC	10.0	.134/3.4	.035	45	ANL0
G47700	KCFW-10.0-35-45-RB-HFANL1-HC	10.0	.134/3.4	.035	45	ANL1
G56279	KCFW-12.0-35-45-RB-HFANL0-HC	12.0	.162/4.1	.035	45	ANL0
G47701	KCFW-12.0-35-45-RB-HFANL1-HC	12.0	.162/4.1	.035	45	ANL1

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Flexor® Ansel Guiding Sheath with Check-Flo Hemostasis Valve

Includes .018" and .035" or .038" compatible dilators

Check-Flo Valve

Order Number	Reference Part Number	Introducer Fr	Sheath ID inch/mm	Accepts Wire Guide Diameter inch	Sheath Length cm	Tip Configuration
G48186	KCFW-4.0-18/35-45-RB-ANL0-HC	4.0	.060/1.5	.018/.035	45	ANL0
G48188	KCFW-4.0-18/35-90-RB-ANL0-HC	4.0	.060/1.5	.018/.035	90	ANL0
G48191	KCFW-4.0-18/35-110-RB-ANL0-HC	4.0	.060/1.5	.018/.035	110	ANL0
G44153	KCFW-5.0-18/38-45-RB-ANL0-HC	5.0	.074/1.9	.018/.038	45	ANL0
G29981	KCFW-5.0-18/38-45-RB-ANL1-HC	5.0	.074/1.9	.018/.038	45	ANL1
G30117	KCFW-5.0-18/38-45-RB-ANL2-HC	5.0	.074/1.9	.018/.038	45	ANL2
G48189	KCFW-5.0-18/38-90-RB-ANL0-HC	5.0	.074/1.9	.018/.038	90	ANL0
G29988	KCFW-5.0-18/38-90-RB-ANL1-HC	5.0	.074/1.9	.018/.038	90	ANL1
G48192	KCFW-5.0-18/38-110-RB-ANL0-HC	5.0	.074/1.9	.018/.038	110	ANL0
G44154	KCFW-6.0-18/38-45-RB-ANL0-HC	6.0	.087/2.2	.038	45	ANL0
G29982	KCFW-6.0-18/38-45-RB-ANL1-HC	6.0	.087/2.2	.038	45	ANL1
G29983	KCFW-6.0-18/38-45-RB-ANL2-HC	6.0	.087/2.2	.038	45	ANL2
G29984	KCFW-6.0-18/38-45-RB-ANL3-HC	6.0	.087/2.2	.038	45	ANL3
G29989	KCFW-6.0-18/38-90-RB-ANL1-HC	6.0	.087/2.2	.038	90	ANL1
G44155	KCFW-7.0-18/38-45-RB-ANL0-HC	7.0	.100/2.5	.038	45	ANL0
G29985	KCFW-7.0-18/38-45-RB-ANL1-HC	7.0	.100/2.5	.038	45	ANL1
G29986	KCFW-7.0-18/38-45-RB-ANL2-HC	7.0	.100/2.5	.038	45	ANL2
G49099	KCFW-8.0-18/38-45-RB-ANL0-HC	8.0	.113/2.9	.038	45	ANL0
G49100	KCFW-8.0-18/38-45-RB-ANL1-HC	8.0	.113/2.9	.038	45	ANL1
G49101	KCFW-9.0-18/38-45-RB-ANL0-HC	9.0	.126/3.2	.038	45	ANL0
G49102	KCFW-9.0-18/38-45-RB-ANL1-HC	9.0	.126/3.2	.038	45	ANL1

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



www.cookmedical.com