


Dauerhafte Symptomlinderung erfordert Ausgewogenheit^{1,6}

Keine Entscheidung
mehr zwischen
Anpassungsfähigkeit
oder Lumenerweiterung¹



Zilver[®] Vena[™]
VENOUS SELF-EXPANDING STENT

Passt sich der Gefäßanatomie an

Das offene Zelldesign sorgt für Flexibilität und minimale Verkürzung, um sich an eine dynamische Umgebung anzupassen.^{1,2,3,4}

Dauerhafte Patency

Das offene Zelldesign und das selbstexpandierende Nitinol-Material sorgen für eine ausreichende Ausdehnung des Lumens ohne unnötigen Kraftaufwand.^{1,2,5}

Kontrollierte Freisetzung

Das Low-Profile-Pin-and-Pull-Platzierungssystem ermöglicht eine präzise und reibungslose Platzierung.^{1,2,3}



Zilver[®] Vena[™]

VENOUS SELF-EXPANDING STENT

- Flexibilität und Knickresistenz der bewährten Zilver[®]-Technologie*
- 2,3 mm (7 French) Schleusenkompatibilität
- Platzierung über 80- und 120-cm-Platzierungssysteme
- Erhöhte Sichtbarkeit durch jeweils vier Goldmarker am proximalen und distalen Ende des Stents**
- Das Pin-and-Pull-Platzierungssystem ermöglicht eine kontrollierte und einfache Platzierung**


* Engineering testing: D00206124 and D00244797

** Engineering testing: D00204576

Offenes Zelldesign



Vier Goldmarker jeweils
am proximalen und
distalen Ende des Stents




Bestellnummer	Referenznummer	Passender Führungsdraht- durchmesser Inch	Stentdurchmesser (mm)	Stentlänge (mm)	Mindestschleuse mm/Fr
80-cm-Over-the-Wire-Platzierungssystem					
G52022	ZVT7-35-80-14-6.0	0,035	14	60	2,3/7
G52023	ZVT7-35-80-14-10.0	0,035	14	100	2,3/7
G52024	ZVT7-35-80-14-14.0	0,035	14	140	2,3/7
G52027	ZVT7-35-80-16-6.0	0,035	16	60	2,3/7
G52026	ZVT7-35-80-16-10.0	0,035	16	100	2,3/7
G52025	ZVT7-35-80-16-14.0	0,035	16	140	2,3/7
120-cm-Over-the-Wire-Platzierungssystem					
G52016	ZVT7-35-120-14-6.0	0,035	14	60	2,3/7
G52018	ZVT7-35-120-14-10.0	0,035	14	100	2,3/7
G52017	ZVT7-35-120-14-14.0	0,035	14	140	2,3/7
G52019	ZVT7-35-120-16-6.0	0,035	16	60	2,3/7
G52021	ZVT7-35-120-16-10.0	0,035	16	100	2,3/7
G52020	ZVT7-35-120-16-14.0	0,035	16	140	2,3/7

Bestimmte Produkte oder Artikelnummern sind eventuell nicht in allen Ländern erhältlich. Einzelheiten erfahren Sie von Ihrem zuständigen Cook-Außendienstmitarbeiter oder beim Cook-Kundendienst.

Lesen Sie die produktspezifischen Informationen zu Risiken in der Gebrauchsanweisung unter cookmedical.eu.

Referenznummer-Schlüssel

ZVT7-35-80-14-6.0



1. Shamimi-Noori SM, Clark TWI. (2018). Venous Stents: Current Status and Future Directions. Tech Vasc Interv Radiol. 2018;21(2):113-116.
2. Zilver® Vena™ View the Instructions for Use for a thorough examination of the procedural instructions, intended use, contraindications, warnings and precautions and potential adverse events. Some products or part numbers may not be available in all markets.
3. Engineering testing: D00206124 and D00204576
4. Engineering testing: D00250996, D00261847, D00250585, and D00244797
5. Engineering testing: D00211273 and D00250794
6. VIVO clinical study. ClinicalTrials.gov Web site. <https://clinicaltrials.gov/ct2/show/results/NCT01970007>
Published October 25, 2013. Updated April 08, 2021. Accessed October 15, 2021.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4