


El alivio duradero de los síntomas requiere equilibrio^{1,6}

Ya no hay que elegir entre adaptabilidad o apertura del lumen¹



Zilver® Vena™
VENOUS SELF-EXPANDING STENT

Se adapta a la anatomía del vaso

El diseño de celda abierta proporciona flexibilidad y acortamiento mínimo para adaptarse a un entorno dinámico.^{1,2,3,4}

Permeabilidad duradera

El diseño de celda abierta y el material autoexpandible de nitinol proporcionan suficiente expansión del lumen sin una fuerza innecesaria.^{1,2,5}

Liberación controlada

El sistema de liberación pin-and-pull de bajo perfil permite una liberación precisa y suave.^{1,2,3}



Zilver[®] Vena[™]

VENOUS SELF-EXPANDING STENT

- Flexibilidad y resistencia al acodamiento de la tecnología probada Zilver[®]*
- Compatibilidad con introductores de 2,3 mm (7 French)
- Despliegue a través de sistemas de liberación de 80 y 120 cm
- Visibilidad mejorada mediante cuatro marcadores de oro en cada extremo proximal y distal del stent**
- El sistema de liberación pin-and-pull permite un despliegue controlado y una fácil colocación**


*Engineering testing: D00206124 and D00244797

**Engineering testing: D00204576

Cuatro marcadores de oro en cada extremo proximal y distal del stent



Diseño de celda abierta




Referencia	Referencia de producto	Acepta diámetro de guía (pulgadas)	Diámetro del stent (mm)	Longitud del stent (mm)	Vaina mínima (mm/Fr)
Sistema de introducción sobre la guía de 80 cm					
G52022	ZVT7-35-80-14-6.0	0,035	14	60	2,3/7
G52023	ZVT7-35-80-14-10.0	0,035	14	100	2,3/7
G52024	ZVT7-35-80-14-14.0	0,035	14	140	2,3/7
G52027	ZVT7-35-80-16-6.0	0,035	16	60	2,3/7
G52026	ZVT7-35-80-16-10.0	0,035	16	100	2,3/7
G52025	ZVT7-35-80-16-14.0	0,035	16	140	2,3/7
Sistema de introducción sobre la guía de 120 cm					
G52016	ZVT7-35-120-14-6.0	0,035	14	60	2,3/7
G52018	ZVT7-35-120-14-10.0	0,035	14	100	2,3/7
G52017	ZVT7-35-120-14-14.0	0,035	14	140	2,3/7
G52019	ZVT7-35-120-16-6.0	0,035	16	60	2,3/7
G52021	ZVT7-35-120-16-10.0	0,035	16	100	2,3/7
G52020	ZVT7-35-120-16-14.0	0,035	16	140	2,3/7

Es posible que algunos productos o números de referencia de productos no estén disponibles en todos los mercados. Póngase en contacto con su representante local de Cook o con el departamento de atención al cliente para obtener más información.

Consulte la información sobre los riesgos del producto en las instrucciones de uso en cookmedical.eu.

Referencia y clave del producto

ZVT7-35-80-14-6.0



1. Shamimi-Noori SM, Clark TWI. (2018). Venous Stents: Current Status and Future Directions. *Tech Vasc Interv Radiol*. 2018;21(2):113-116.
2. Zilver® Vena™ View the Instructions for Use for a thorough examination of the procedural instructions, intended use, contraindications, warnings and precautions and potential adverse events. Some products or part numbers may not be available in all markets.
3. Engineering testing: D00206124 and D00204576
4. Engineering testing: D00250996, D00261847, D00250585, and D00244797
5. Engineering testing: D00211273 and D00250794
6. VIVO clinical study. ClinicalTrials.gov Web site. <https://clinicaltrials.gov/ct2/show/results/NCT01970007>. Published October 25, 2013. Updated April 08, 2021. Accessed October 15, 2021.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4