



Accessing opportunity^{1,2}

Designed to support device technology used to access and cross vascular obstructions^{1,2,3}



Four configurations and two lengths available



TriForce[®] PERIPHERAL CROSSING SET

Disease and anatomy specific


Designed to support device technology used to access and cross vascular obstructions.^{1,2,3}

Responsive torque

Bringing together the inherent columnar strength of the Flexor[®] sheath and the stainless steel braid of the CXI[®] crossing catheter, allowing continuous advancement.

Enhanced pushability and trackability

Hydrophilically coated and tungsten loaded tip to enhance crossing ability.



TriForce®

PERIPHERAL CROSSING SET

- The TriForce Peripheral Crossing Set includes a 5.0 Fr Flexor outer sheath, a 4.0 Fr CXI Support Catheter, and a Peel-Away® sheath.
- The Flexor outer sheath tapers to a 4.0 Fr inner diameter at the tip to ensure a smooth transition to the CXI Support Catheter.
- The Flexor outer sheath features radiopaque tip and hydrophilic coating.
- TriForce accepts a 0.035 inch diameter wire guide.

Order Number	Reference Part Number	Catheter Fr/Length cm	Sheath Length cm	Sheath Tip Fr	Accepts Wire Guide Diameter Inch
Set with Straight-Tip Sheath and Straight-Tip Catheter					
G56412	KCXS-5.0-35-65-RB-0/0-HC	4.0/65	55	4.0	0.035
G56416	KCXS-5.0-35-100-RB-0/0-HC	4.0/100	90	4.0	0.035
Set with Straight-Tip Sheath and Curved-Tip Catheter					
G56413	KCXS-5.0-35-65-RB-0/DAV-HC	4.0/65	55	4.0	0.035
G56417	KCXS-5.0-35-100-RB-0/DAV-HC	4.0/100	90	4.0	0.035
Set with Curved-Tip Sheath and Straight-Tip Catheter					
G56414	KCXS-5.0-35-65-RB-MPB/0-HC	4.0/65	55	4.0	0.035
G56418	KCXS-5.0-35-100-RB-MPB/0-HC	4.0/100	90	4.0	0.035
Set with Curved-Tip Sheath and Curved-Tip Catheter					
G56415	KCXS-5.0-35-65-RB-MPB/DAV-HC	4.0/65	55	4.0	0.035
G56419	KCXS-5.0-35-100-RB-MPB/DAV-HC	4.0/100	90	4.0	0.035

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Please see product risk information in the IFU at cookmedical.eu.

1. Hofmann L. What's new and what's needed in chronic venous occlusion therapy. Endovascular Today Web site. <https://evtoday.com/articles/2015-oct/whats-new-and-whats-needed-in-chronic-venous-occlusion-therapy>. Published October 2015. Accessed December 20, 2021.
2. Case study: Erbacher G. Case study: Crossing a difficult venous occlusion with the TriForce® Peripheral Crossing Set, Oklahoma State University Medical Center, PI-D64062-EN. https://www.cookmedical.com/data/resources/PI-D52143-EN-F_M3_1576703674956.pdf
3. Lu T, Schwein A, Bechara C, Bismuth J. Improved percutaneous Central venous recanalization: rates with novel telescoping catheter system. JVS 2016;63(6) (suppl):132S-133S. [https://www.jvascsurg.org/article/S0741-5214\(16\)00539-5/fulltext](https://www.jvascsurg.org/article/S0741-5214(16)00539-5/fulltext)

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4