

# Ensure an effective filter retrieval<sup>1,2,3</sup>



## Günther Tulip®

VENA CAVA FILTER RETRIEVAL SET FOR JUGULAR VEIN APPROACH

### Successful retrieval

The successful retrieval rate was 90.2%<sup>2</sup> when removing Günther Tulip Vena Cava Filters, and the successful retrieval rate was 96.6%<sup>3</sup> when removing Celect® Vena Cava Filters.

### Designed for ease-of-use

The single loop design was created and tested for the retrieval of Cook filters. The snare loop size fits the vena cava to facilitate the capture of the filter hook.

### Solution in a complete set

The retrieval set contains snare, coaxial sheath set, access needle, and wire guide.



## Günther Tulip®

VENA CAVA FILTER RETRIEVAL SET FOR JUGULAR VEIN APPROACH

Used for the retrieval of implanted Günther Tulip and Cook Celect Platinum Vena Cava Filters in patients who no longer require a filter. Retrieval of the filter can be performed only by a jugular approach.

Order Number	Reference Part Number	Retrieval Loop System Catheter Diameter Fr	Retrieval Loop System Length cm	Sheath Diameter Fr	Sheath Working Length cm
G13287	GTRS-200-RB	6.3	80	11.0	60

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

1. Hoppe H, Nutting CW, Smouse HR, et al. Günther Tulip filter retrievability multicenter study including CT follow-up: final report. *J Vasc Interv Radiol.* 2006;17(6):1017-1023.
2. Smouse HB, Rosenthal D, Thuong VH, et al. Long-term retrieval success rate profile for the Günther Tulip vena cava filter. *J Vasc Interv Radiol.* 2009;20(7):871-877.
3. Lyon SM, Riojas GE, Uberoi R, et al. Short- and long-term retrievability of the Celect vena cava filter: results from a multi-institutional registry. *J Vasc Interv Radiol.* 2009;20(11):1441-1448.

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 800 7615, [is.orders@cookmedical.com](mailto:is.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)

EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)

#### Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)

#### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI-ESC-IR-OHNS-PI-RH-SUR-A4