


Improved crossability¹



Starcast™ shoulders

Pliiform™ balloon texture

Advance® 35LP

LOW-PROFILE PTA BALLOON DILATATION CATHETER

Designed for reinsertion


Starcast balloon shoulders facilitate sheath reinsertion.

Pliant surface texture

Patented balloon surface texture helps you cross lesions.¹

Trackability

Flexible tip enhances your ability to navigate tortuous vasculature.²




**D-shaped lumen design
for inflating and deflating balloon**

1. Burton D, Lysgaard T, Aggerholm S, et al. Balloon catheter and method for making same. United States patent Pub No. US20100036314A1, published 11 February 2010.

2. Animal evaluation 07-19-2012 (Purdue University / MED Institute): Advance 35LP line ext. PTA catheter. DV1-072512AC (2012).

ORDERING INFORMATION

Reference Part Number key



Order Number	Reference Part Number	Wire Guide Diameter inch	Balloon Inflated Diameter mm	Balloon Length cm	Recommended Introducer Fr	Nominal pressure atm	Rated Burst Pressure atm
0.035 inch Diameter Wire Guide – 80 cm Long Shaft							
G52325	PTA5-35-80-3-2.0	0.035	3	2	5.0	10	15
G52326	PTA5-35-80-3-4.0	0.035	3	4	5.0	10	15
G52327	PTA5-35-80-3-8.0	0.035	3	8	5.0	10	15
G52328	PTA5-35-80-3-10.0	0.035	3	10	5.0	10	15
G52330	PTA5-35-80-3-14.0	0.035	3	14	5.0	10	15
G34350	PTA5-35-80-3-20.0	0.035	3	20	5.0	10	15
G52331	PTA5-35-80-4-2.0	0.035	4	2	5.0	10	15
G52332	PTA5-35-80-4-4.0	0.035	4	4	5.0	10	15
G34355	PTA5-35-80-4-6.0	0.035	4	6	5.0	10	15
G52333	PTA5-35-80-4-8.0	0.035	4	8	5.0	10	15
G52334	PTA5-35-80-4-10.0	0.035	4	10	5.0	10	15
G52336	PTA5-35-80-4-14.0	0.035	4	14	5.0	10	15
G34354	PTA5-35-80-4-20.0	0.035	4	20	5.0	10	15
G52337	PTA5-35-80-5-2.0	0.035	5	2	5.0	10	15
G52338	PTA5-35-80-5-4.0	0.035	5	4	5.0	10	15
G52339	PTA5-35-80-5-6.0	0.035	5	6	5.0	10	15
G52340	PTA5-35-80-5-8.0	0.035	5	8	5.0	10	15
G52341	PTA5-35-80-5-10.0	0.035	5	10	5.0	10	15
G52343	PTA5-35-80-5-14.0	0.035	5	14	5.0	10	15
G34359	PTA5-35-80-5-20.0	0.035	5	20	5.0	10	15
G52344	PTA5-35-80-6-2.0	0.035	6	2	5.0	10	15
G52345	PTA5-35-80-6-4.0	0.035	6	4	5.0	10	15
G52346	PTA5-35-80-6-6.0	0.035	6	6	5.0	10	15
G52347	PTA5-35-80-6-8.0	0.035	6	8	5.0	10	15
G52348	PTA5-35-80-6-10.0	0.035	6	10	5.0	10	15
G52350	PTA5-35-80-6-14.0	0.035	6	14	5.0	10	15
G34363	PTA5-35-80-6-20.0	0.035	6	20	5.0	10	15
G52351	PTA5-35-80-7-2.0	0.035	7	2	5.0	8	12
G52352	PTA5-35-80-7-4.0	0.035	7	4	5.0	8	12
G52353	PTA5-35-80-7-6.0	0.035	7	6	5.0	8	12
G52354	PTA5-35-80-7-8.0	0.035	7	8	5.0	8	12
G52355	PTA5-35-80-7-10.0	0.035	7	10	6.0	8	12
G52357	PTA5-35-80-7-14.0	0.035	7	14	6.0	8	12
G34367	PTA5-35-80-7-20.0	0.035	7	20	6.0	8	12

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Continued on next page

ORDERING INFORMATION

Order Number	Reference Part Number	Wire Guide Diameter inch	Balloon Inflated Diameter mm	Balloon Length cm	Recommended Introducer Fr	Nominal pressure atm	Rated Burst Pressure atm
0.035 inch Diameter Wire Guide – 80 cm Long Shaft							
G52358	PTA5-35-80-8-2.0	0.035	8	2	5.0	8	11
G52359	PTA5-35-80-8-4.0	0.035	8	4	5.0	8	11
G52360	PTA5-35-80-8-6.0	0.035	8	6	5.0	8	11
G52361	PTA5-35-80-8-8.0	0.035	8	8	5.0	8	11
G52362	PTA5-35-80-8-10.0	0.035	8	10	6.0	8	11
G52364	PTA5-35-80-8-14.0	0.035	8	14	6.0	8	11
G34371	PTA5-35-80-8-20.0	0.035	8	20	6.0	8	11
G34372	PTA5-35-80-9-2.0	0.035	9	2	6.0	8	11
G52365	PTA5-35-80-9-3.0	0.035	9	3	6.0	8	11
G52366	PTA5-35-80-9-4.0	0.035	9	4	6.0	8	11
G52367	PTA5-35-80-9-6.0	0.035	9	6	6.0	8	11
G52368	PTA5-35-80-9-8.0	0.035	9	8	6.0	8	11
G52369	PTA5-35-80-9-10.0	0.035	9	10	6.0	8	11
G52371	PTA5-35-80-9-14.0	0.035	9	14	6.0	8	11
G34373	PTA5-35-80-10-2.0	0.035	10	2	6.0	8	11
G52372	PTA5-35-80-10-3.0	0.035	10	3	6.0	8	11
G52373	PTA5-35-80-10-4.0	0.035	10	4	6.0	8	11
G52374	PTA5-35-80-10-6.0	0.035	10	6	6.0	8	11
G52375	PTA5-35-80-10-8.0	0.035	10	8	6.0	8	11
G52376	PTA5-35-80-10-10.0	0.035	10	10	7.0	8	11
G52378	PTA5-35-80-10-14.0	0.035	10	14	7.0	8	11
G34374	PTA5-35-80-12-2.0	0.035	12	2	7.0	5	8
G52379	PTA5-35-80-12-4.0	0.035	12	4	7.0	5	8
G52380	PTA5-35-80-12-8.0	0.035	12	8	7.0	5	8
G52381	PTA5-35-80-12-10.0	0.035	12	10	7.0	5	8
G52383	PTA5-35-80-12-14.0	0.035	12	14	7.0	5	8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Order Number	Reference Part Number	Accepts Wire Guide Diameter inch	Balloon Inflated Diameter mm	Balloon Length cm	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
0.035 inch Diameter Wire Guide – 135 cm Long Shaft							
G52245	PTA5-35-135-3-2.0	0.035	3	2	5.0	10	15
G52246	PTA5-35-135-3-4.0	0.035	3	4	5.0	10	15
G34375	PTA5-35-135-3-6.0	0.035	3	6	5.0	10	15
G52247	PTA5-35-135-3-8.0	0.035	3	8	5.0	10	15
G52248	PTA5-35-135-3-10.0	0.035	3	10	5.0	10	15
G52249	PTA5-35-135-3-12.0	0.035	3	12	5.0	10	15
G52250	PTA5-35-135-3-14.0	0.035	3	14	5.0	10	15
G34377	PTA5-35-135-3-20.0	0.035	3	20	5.0	10	15
G52251	PTA5-35-135-4-2.0	0.035	4	2	5.0	10	15
G52252	PTA5-35-135-4-4.0	0.035	4	4	5.0	10	15
G35528	PTA5-35-135-4-6.0	0.035	4	6	5.0	10	15
G52253	PTA5-35-135-4-8.0	0.035	4	8	5.0	10	15
G52254	PTA5-35-135-4-10.0	0.035	4	10	5.0	10	15
G52255	PTA5-35-135-4-12.0	0.035	4	12	5.0	10	15
G52256	PTA5-35-135-4-14.0	0.035	4	14	5.0	10	15
G35527	PTA5-35-135-4-20.0	0.035	4	20	5.0	10	15
G52257	PTA5-35-135-5-2.0	0.035	5	2	5.0	10	15
G52258	PTA5-35-135-5-4.0	0.035	5	4	5.0	10	15
G52259	PTA5-35-135-5-6.0	0.035	5	6	5.0	10	15
G52260	PTA5-35-135-5-8.0	0.035	5	8	5.0	10	15
G52261	PTA5-35-135-5-10.0	0.035	5	10	5.0	10	15
G52262	PTA5-35-135-5-12.0	0.035	5	12	5.0	10	15
G52263	PTA5-35-135-5-14.0	0.035	5	14	5.0	10	15
G35531	PTA5-35-135-5-17.0	0.035	5	17	5.0	10	15
G35532	PTA5-35-135-5-20.0	0.035	5	20	5.0	10	15
G52264	PTA5-35-135-6-2.0	0.035	6	2	5.0	10	15
G52265	PTA5-35-135-6-4.0	0.035	6	4	5.0	10	15
G52266	PTA5-35-135-6-6.0	0.035	6	6	5.0	10	15
G52267	PTA5-35-135-6-8.0	0.035	6	8	5.0	10	15
G52268	PTA5-35-135-6-10.0	0.035	6	10	5.0	10	15
G52269	PTA5-35-135-6-12.0	0.035	6	12	5.0	10	15
G52270	PTA5-35-135-6-14.0	0.035	6	14	5.0	10	15
G35535	PTA5-35-135-6-17.0	0.035	6	17	5.0	10	15
G35536	PTA5-35-135-6-20.0	0.035	6	20	5.0	10	15
G52271	PTA5-35-135-7-2.0	0.035	7	2	5.0	8	12
G52272	PTA5-35-135-7-4.0	0.035	7	4	5.0	8	12
G52273	PTA5-35-135-7-6.0	0.035	7	6	5.0	8	12
G52274	PTA5-35-135-7-8.0	0.035	7	8	5.0	8	12
G52275	PTA5-35-135-7-10.0	0.035	7	10	6.0	8	12
G52276	PTA5-35-135-7-12.0	0.035	7	12	6.0	8	12
G52277	PTA5-35-135-7-14.0	0.035	7	14	6.0	8	12
G35540	PTA5-35-135-7-20.0	0.035	7	14	6.0	8	12

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

Continued on next page

ORDERING INFORMATION











Order Number	Reference Part Number	Accepts Wire Guide Diameter inch	Balloon Inflated Diameter mm	Balloon Length cm	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
0.035 inch Diameter Wire Guide – 135 cm Long Shaft							
G52278	PTA5-35-135-8-2.0	0.035	8	2	5.0	8	11
G52279	PTA5-35-135-8-4.0	0.035	8	4	5.0	8	11
G52280	PTA5-35-135-8-6.0	0.035	8	6	5.0	8	11
G52281	PTA5-35-135-8-8.0	0.035	8	8	5.0	8	11
G52282	PTA5-35-135-8-10.0	0.035	8	10	6.0	8	11
G52283	PTA5-35-135-8-12.0	0.035	8	12	6.0	8	11
G52284	PTA5-35-135-8-14.0	0.035	8	14	6.0	8	11
G35544	PTA5-35-135-8-20.0	0.035	8	20	6.0	8	11
G35545	PTA5-35-135-9-2.0	0.035	9	2	6.0	8	11
G52285	PTA5-35-135-9-3.0	0.035	9	3	6.0	8	11
G52286	PTA5-35-135-9-4.0	0.035	9	4	6.0	8	11
G52287	PTA5-35-135-9-6.0	0.035	9	6	6.0	8	11
G52288	PTA5-35-135-9-8.0	0.035	9	8	6.0	8	11
G52289	PTA5-35-135-9-10.0	0.035	9	10	6.0	8	11
G52291	PTA5-35-135-9-14.0	0.035	9	14	6.0	8	11
G35546	PTA5-35-135-10-2.0	0.035	10	2	6.0	8	11
G52292	PTA5-35-135-10-3.0	0.035	10	3	6.0	8	11
G52293	PTA5-35-135-10-4.0	0.035	10	4	6.0	8	11
G52294	PTA5-35-135-10-6.0	0.035	10	6	6.0	8	11
G52295	PTA5-35-135-10-8.0	0.035	10	8	6.0	8	11
G52296	PTA5-35-135-10-10.0	0.035	10	10	7.0	8	11
G52298	PTA5-35-135-10-14.0	0.035	10	14	7.0	8	11
G35547	PTA5-35-135-12-2.0	0.035	12	2	7.0	5	8
G52299	PTA5-35-135-12-4.0	0.035	12	2	7.0	5	8
G52300	PTA5-35-135-12-8.0	0.035	12	8	7.0	5	8
G52301	PTA5-35-135-12-10.0	0.035	12	10	7.0	5	8
G52303	PTA5-35-135-12-14.0	0.035	12	14	7.0	5	8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.

		Balloon Diameter (mm)								
		3	4	5	6	7	8	9	10	12
Pressure (bar)(atm)	4	2.64	3.49	4.28	5.26	6.31	7.17	8.08	9.11	11.62
	5	2.69	3.58	4.40	5.39	6.52	7.42	8.37	9.40	12.00
	6	2.75	3.67	4.54	5.55	6.72	7.66	8.62	9.65	12.32
	7	2.82	3.77	4.67	5.68	6.88	7.85	8.84	9.84	12.58
	8	2.87	3.86	4.79	5.80	7.00	8.00	9.00	10.00	12.84
	9	2.95	3.93	4.90	5.90	7.12	8.15	9.15	10.16	
	10	3.00	4.00	5.00	6.00	7.23	8.27	9.33	10.29	
	11	3.05	4.06	5.09	6.09	7.34	8.41	9.49	10.47	
	12	3.10	4.12	5.19	6.17	7.46				
	13	3.15	4.18	5.28	6.26					
	14	3.20	4.24	5.37	6.34					
	15	3.25	4.31	5.47	6.44					

■ Nominal Balloon Pressure
■ Rated Burst Pressure (RBP)

Providing you procedural support where it matters.

ACCESS	TARGET	CROSS	TREAT
<p data-bbox="239 504 383 548">Micropuncture[®] Introducer Set</p> 	<p data-bbox="550 504 718 548">Amplatz Support Wire Guide</p> 	<p data-bbox="893 504 1021 548">CXI[®] Support Catheter</p> 	<p data-bbox="1133 504 1428 548">Advance[®] 35LP Low-Profile PTA Balloon Dilatation Catheter</p> 
<p data-bbox="215 896 406 940">Check-Flo Performer[®] Introducer</p> 	<p data-bbox="526 896 742 940">Roadrunner[®] UniGlide[®] Hydrophilic Wire Guide</p> 		<p data-bbox="1165 896 1396 940">Zilver[®] PTX[®] Drug-Eluting Peripheral Stent</p> 
<p data-bbox="255 1299 367 1344">Bentson Wire Guide</p> 	<p data-bbox="534 1321 734 1344">Beacon[®] Tip Catheter</p> 		
	<p data-bbox="526 1691 742 1736">Flexor[®] Introducer / Flexor[®] Guiding Sheath</p> 		

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4