



Advance Serenity 14

HYDROPHILIC PTA BALLOON DILATATION CATHETER

Pressure (atm)	Balloon Diameter (mm)				
	2	2.5	3	3.5	4
2	1.75	2.15	2.50	2.95	2.08
4	1.81	2.40	2.86	3.39	3.83
6	1.91	2.52	3.00	3.58	4.02
8	2.01	2.64	3.12	3.73	4.18
10	2.09	2.73	3.21	3.85	4.31
12	2.15	2.80	3.28	3.94	4.40
14	2.19	2.85	3.34	4.01	4.48
16	2.23	2.90	3.39	4.07	4.54

Nominal pressure

Rated burst pressure

Advance Serenity 18

HYDROPHILIC PTA BALLOON DILATATION CATHETER

Pressure (atm)	Balloon Diameter (mm)											
	2	2.5	3	3.5	4	5	6	7	8 (2-4 cm length)	8 (6-8 cm length)	9	10
8	1.91	2.41	2.92	3.51	4.00	4.85	5.80	6.83	7.83	7.95	8.87	9.95
10	1.98	2.50	3.04	3.63	4.14	5.01	5.99	7.01	8.02	8.14	9.09	10.21
12	2.04	2.57	3.13	3.73	4.25	5.12	6.13	7.14	8.17	8.30		
14	2.09	2.63	3.20	3.80	4.33	5.22	6.24	7.25				
16	2.13	2.67	3.25	3.86	4.41	5.33	6.35					