



Accessing opportunity

Access and cross vascular obstructions with just one set^{1,2,3}



Four configurations and two lengths available



TriForce[®] PERIPHERAL CROSSING SET

Designed for responsive torque

TriForce brings together the inherent columnar strength of the Flexor[®] sheath and the stainless steel braid of the CXI[®] crossing catheter, allowing continuous advancement.⁴

Visibility

Distal ends of the catheter and sheath are designed with increased radiopacity for visualization.⁴

Optimized trackability

Designed coaxial pairing of CXI catheter with Flexor sheath facilitates advancement. Hydrophilically coated and tungsten loaded tip is designed to enhance crossing ability.⁴




1. Hofmann L. What's new and what's needed in chronic venous occlusion therapy. Endovascular Today Web site.

<https://evtoday.com/articles/2015-oct/whats-new-and-whats-needed-in-chronic-venous-occlusion-therapy>. Published in October 2015. Accessed on December 20, 2021.

2. Case study: Crossing a difficult venous occlusion with the TriForce[®] Peripheral Crossing Set, George Erbächer, DO, Oklahoma State University Medical Center, PI-D52143-EN.

3. Lu T, Schwein A, Bechara C, et al. Improved percutaneous central venous recanalization rates with novel telescoping catheter system. *J Vasc Surg.* 2016;63(6)(suppl):1325-1335.

4. TriForce[®] Peripheral Crossing Set. View the Instructions for Use for a thorough examination of the procedural instructions, intended use, contraindications, warnings and precautions, and potential adverse events.



TriForce®

PERIPHERAL CROSSING SET

- TriForce Peripheral Crossing Set includes an outer 5.0 Fr Flexor sheath, a 4.0 Fr CXI Support Catheter, and a Peel-Away® sheath.
- Flexor outer sheath tapers to a 4.0 Fr inner diameter at the tip to ensure a smooth transition to the CXI Support Catheter.
- Flexor outer sheath features radiopaque tip and hydrophilic coating.
- TriForce accepts a 0.035 inch diameter wire guide.

Order Number	Reference Part Number	Catheter Fr/Length cm	Sheath ID/Length Fr/cm	Sheath Tip Fr	Accepts Wire Guide Diameter Inch
Set with Straight Tip Sheath and Straight Tip Catheter					
G56412	KCXS-5.0-35-65-RB-0/0-HC	4.0/65	5.0/55	4.0	0.035
G56416	KCXS-5.0-35-100-RB-0/0-HC	4.0/100	5.0/90	4.0	0.035
Set with Straight Tip Sheath and Curved Tip Catheter					
G56413	KCXS-5.0-35-65-RB-0/DAV-HC	4.0/65	5.0/55	4.0	0.035
G56417	KCXS-5.0-35-100-RB-0/DAV-HC	4.0/100	5.0/90	4.0	0.035
Set with Curved Tip Sheath and Straight Tip Catheter					
G56414	KCXS-5.0-35-65-RB-MPB/0-HC	4.0/65	5.0/55	4.0	0.035
G56418	KCXS-5.0-35-100-RB-MPB/0-HC	4.0/100	5.0/90	4.0	0.035
Set with Curved Tip Sheath and Curved Tip Catheter					
G56415	KCXS-5.0-35-65-RB-MPB/DAV-HC	4.0/65	5.0/55	4.0	0.035
G56419	KCXS-5.0-35-100-RB-MPB/DAV-HC	4.0/100	5.0/90	4.0	0.035

Set includes 5.0 Fr Flexor outer sheath, 4.0 Fr CXI Support Catheter and Peel-Away sheath.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11