


# Reliable PE prevention with Celect Platinum<sup>®1,2</sup>

1.55% recurrent PE rate<sup>1,2</sup>

Effective at preventing PE while also being highly retrievable.<sup>1</sup>



Platinum marker



## Celect Platinum<sup>®</sup> VENA CAVA FILTER

### Clinically successful design

The four radiopaque platinum markers can verify placement accuracy.<sup>3</sup> Anchors and platinum markers on the filter's feet allow secure caval fixation.<sup>4</sup>

### Accurate, controlled deployment

The NavAlign<sup>™</sup> delivery system helps to deploy the filter precisely, easily providing control with accurate deployment.<sup>5</sup>

### Constructed of unique material

Cook filters are constructed of a durable Conichrome<sup>®</sup> alloy.<sup>1,6</sup>



# Celect Platinum®

VENA CAVA FILTER

The Cook Celect Platinum Filter implant is indicated for prevention of symptomatic pulmonary embolism (PE) by capturing blood clots travelling in the infrarenal IVC in the following clinical situations:

- Patients with evidence of PE or deep vein thrombosis (DVT) in the IVC or lower extremity veins and either an absolute or relative contraindication to anticoagulation, a failure of anticoagulation, or high risk from complication from anticoagulation
- Patients in the perioperative setting with a history of PE or DVT for whom anticoagulation must be interrupted
- Patients at high risk of DVT or PE but considered unable to receive anticoagulation due to a high risk of complication

The Cook Celect Platinum Filter implant may be retrieved if clinically indicated; please refer to the “Optional Filter Retrieval” section of the Instructions for Use for more information.

Order Number	Reference Part Number	Sheath Diameter Fr	Max. Indicated Cava Diameter for Filter Placement mm	Filter Length mm	Sheath Working Length cm
<b>Celect Platinum® Vena Cava Filter Set for Femoral Vein Approach</b>					
G34501	IGTCFS-65-2-FEM-CELECT-PT	7	30	49	65
<b>Celect Platinum® Vena Cava Filter Set for Jugular Vein Approach</b>					
G34310	IGTCFS-65-2-JUG-CELECT-PT	7	30	49	65
<b>Celect Platinum® Vena Cava Filter Set for Femoral and Jugular Vein Approach</b>					
G34504	IGTCFS-65-2-UNI-CELECT-PT	7	30	49	65

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Please see product risk information in the IFU at [cookmedical.eu](http://cookmedical.eu).

1. Lyon SM, Riojas GE, Uberoi R, et al. Short- and long-term retrievability of the Celect vena cava filter: results from a multi-institutional registry. *J Vasc Interv Radiol.* 2009;20(11):1441-1448.
2. Prospective study of the Cook Celect Filter, including permanent and retrievable use (Study Protocol #05-507) final study report; September 2009.
3. Demonstrated in simulated use tests: D00302649 and D00006045 (TR 4524).
4. Demonstrated in testing: D00006050 (TR 3789) and D00006018 (TR-3787).
5. Celect Platinum® Vena Cava Filter: View the Instructions for Use for a thorough examination of the procedural instructions, intended use, contraindications, warnings and precautions, and potential adverse events.
6. Lee BE, Van Allan RJ, Friedman ML, et al. Complications and retrieval characteristics of Celect Platinum inferior vena cava filters. *J Vasc Surg Venous Lymphat Disord.* 2018;6(2):163-172.

## Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Iceland: +354 800 7615, [is.orders@cookmedical.com](mailto:is.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)  
 EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)  
**Americas:**  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)  
**Australia:**  
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI-ESC-IR-OHNS-PI-RH-SUR-A4