

Rösch-Thurmond


FALLOPIAN TUBE CATHETERIZATION SET

Intended to be used through a uterine access device for fluoroscopic selective catheterization of the proximal fallopian tubes and injection of contrast medium in the evaluation of tubal patency.

- Specifically indicated and designed for the proximal fallopian tubes
- Includes all components for tubal catheterization
- Designed to be placed through uterine access devices that have 6 Fr or larger working channels



The Rösch-Thurmond Fallopian Tube Catheterization Set contains these components:



Order Number	Reference Part Number	Torque Control Catheter Fr	Injection Catheter Fr	Cope Mandril Wire Guide Diameter inch/mm	Safe-T-J Wire Guide Diameter inch/mm
Rösch-Thurmond Fallopian Tube Catheterization Set					
G04324	FTC-900	5.5	3.0	0.015/0.38	0.035/0.89
Accessory Order Number	Reference Part Number		Tip Fr	Length cm	Silicone Balloon Volume mL
Intrauterine Access Balloon Catheter					
G30397	J-IAB-090023		9.0	23	2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Customer Service

EU Website: cookmedical.eu

EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001702, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com Italv: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Norway: +47 23162968, no.orders@cookmedical.com

Iceland: +354 800 7615, is.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com EDI: cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



