

Manage postpartum hemorrhage.




Bakri®

POSTPARTUM BALLOON WITH RAPID INSTILLATION COMPONENTS

A simple solution for postpartum hemorrhage (PPH)

Sequential treatment considerations for conservative management of uterotonic-resistant PPH^{1,2,3}



1. California Maternal Quality Care Collaborative. Obstetric hemorrhage emergency management plan: table chart format. Version 2.0. Published 2014. Accessed August 31, 2016.

2. Committee on Practice Bulletins-Obstetrics. Practice bulletin no. 183: postpartum hemorrhage. *Obstet Gynecol.* 2017;130(4):e168-e186.

3. <https://physicians.utah.edu/echo/pdfs/pregnancy-care-didactics/2016-09.09-postpartum-hemorrhage.pdf> (accessed 29th October 2020)


The Bakri Postpartum Balloon with Rapid Instillation Components is used for temporary control or reduction of postpartum uterine bleeding when conservative management of uterine bleeding is warranted.



When inflated with sterile liquid, the Bakri Postpartum Balloon with Rapid Instillation Components applies pressure to the uterine walls that may provide management of postpartum hemorrhage.

The balloon and catheter are constructed of 100% silicone.

- Can be used after vaginal or cesarean delivery
- Includes rapid instillation components to facilitate inflation of the balloon
- May be used with B-Lynch compression sutures if clinically warranted






Bakri®

POSTPARTUM BALLOON WITH RAPID INSTILLATION COMPONENTS

This device is intended to provide temporary control or reduction of postpartum uterine bleeding when conservative management is warranted. The set includes a syringe, 180 cm long extension tubing with bag spike, and dual check valve. Please refer to the Instructions for Use for complete information on product usage, indications, and contraindications.



Order Number	Reference Part Number	Catheter Fr	Length cm	Balloon Volume mL
G24237	J-SOSR-100500	24	54	500

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Scan this QR code to see other resources.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4