


Inject with confidence.



Williams

CYSTOSCOPIC INJECTION NEEDLE



Williams

CYSTOSCOPIC INJECTION NEEDLE

The Williams Cystoscopic Injection Needle is used for cystoscopic injection into the urethra, bladder neck, and bladder wall.



L'ago per iniezione cistoscopica Williams è utilizzato per l'iniezione cistoscopica nell'uretra, nel collo vescicale e nella parete della vescica.



Игла Williams для цистоскопической инъекции применяется для выполнения цистоскопической инъекции в мочеиспускательный канал, шейку и стенку мочевого пузыря.



Die Williams Zystoskop-Injektionskanüle wird zur Injektion in die Urethra, den Blasenhal und die Blasenwand mittels Zystoskop verwendet.



La aguja de inyección cistoscópica Williams se utiliza para la inyección cistoscópica en la uretra, el cuello vesical y la pared vesical.



Cystoskopowa igła Williams do wstrzyknięć służy do wstrzykiwania do cewki moczowej, szyi i ściany pęcherza moczowego pod kontrolą cystoskopu.

Order Number	Reference Part Number	Fr	Gage	Length cm	Tip Length mm	Comments
G15296	090001-S10	3.7	23	35	8	thin-wall cannula
G16684	090001-S34	3.7	23	35	8	extra-thin-wall cannula
G14220	090001	5	23	35	8	-
G16112	090001-S26	5	23	45	8	-
G15276	090001-S9	5	25	35	8	-

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



CC, URO-A4