



Look at resection from a **variety of angles.**

Bipolar Resectoscopic Electrodes




Bipolar Transurethral Cutting Loop




| Order Number | Reference Part Number | Product Description | Wire Angle | Sheath Fr |
|--------------|-----------------------|--------------------------|------------|-----------|
| G59108 | BTCL-2412-LG | 12° Cutting Loop, Large | 12° | 24 |
| G59109 | BTCL-2412-MD | 12° Cutting Loop, Medium | 12° | 24 |
| G59110 | BTCL-2430-LG | 30° Cutting Loop, Large | 30° | 24 |
| G59111 | BTCL-2430-MD | 30° Cutting Loop, Medium | 30° | 24 |

Bipolar Transurethral Bladder Loop




| Order Number | Reference Part Number | Product Description | Wire Angle | Sheath Fr |
|--------------|-----------------------|---------------------|------------|-----------|
| G59112 | BTBL-2400 | Bladder Loop | 136° | 24 |

Bipolar Transurethral Needle Electrode



| Order Number | Reference Part Number | Product Description | Wire Angle | Sheath Fr |
|--------------|-----------------------|---------------------|------------|-----------|
| G59113 | BTNE-2400 | Needle Electrode | 135° | 24 |

Bipolar Transurethral Plasma Disc®



| Order Number | Reference Part Number | Product Description | Disc Diameter inch | Sheath Fr |
|--------------|-----------------------|--------------------------|--------------------|-----------|
| G59114 | BTPD-2400 | Plasma Disc Vaporization | 0.125 | 24 |

Sold in boxes of 6. Designed to fit the Olympus bipolar WA22366A active and the WA22367A passive working elements. Use with Olympus PK Superpulse, Olympus UES-40, and Olympus ESG-400 generator systems. Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Plasma Disc is a registered trademark of Omnittech Systems Inc.

Customer Service

EU Website: cookmedical.eu

EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001702, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Iceland: +354 8007615, IS.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222

E-mail: cau.custserv@cookmedical.com