


Für die Harnableitung



Bander Ureteral Diversion Stent Set





Bänder Ureteral Diversion Stent Set

Zur intraoperativen Platzierung eines Ureterstents und zur Drainage aus der Niere in das externe Stoma während der Anlage eines uretero-ilealen Leitungswegs und zur kontinenten Harnableitung.

Der Stent darf nicht länger als 90 Tage im Körper verweilen.

Das Set enthält zwei röntgendichte Bänder Stents, zwei Katheterhalter, einen Katheteradapter und einen Führungsdraht (0,038 Inch).



Bestellnummer	Referenznummer	Fr	Länge cm
Set			
G16997	025226	6,0	75
G16996	025227	7,2	75
G16995	025228	8,4	75
Offenes Ende			
G17186	025226-OE	6,0	75
G18070	025227-OE	7,2	75
G17185	025228-OE	8,4	75
Rechter Ureter			
G14779	025707	7,2	75
G16187	025708	8,4	75
Rechter Ureter - offenes Ende			
G16818	025706-S1	6,0	75
G14822	025707-S1	7,2	70
Linker Ureter			
G14780	025807	7,2	75
Linker Ureter - offenes Ende			
G16817	025806-S1	6,0	75
G14823	025807-S1	7,2	70

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4

