



Create a diversion



Bander Ureteral Diversion Stent Set




Bander Ureteral Diversion Stent Set

Used for intraoperative placement to stent the ureter and provide drainage from the kidney to the external stoma during ureteroileal conduit construction and continent urinary diversions.

The stent must not remain indwelling more than 90 days.

The set includes two radiopaque Bander stents, two catheter retainers, a catheter adapter, and a 0.038 inch wire guide.



Order Number	Reference Part Number	Fr	Length cm
Set			
G16997	025226	6.0	75
G16996	025227	7.2	75
G16995	025228	8.4	75
Open End			
G17186	025226-OE	6.0	75
G18070	025227-OE	7.2	75
G17185	025228-OE	8.4	75
Right Ureter			
G14779	025707	7.2	75
G16187	025708	8.4	75
Right Ureter-Open End			
G16818	025706-S1	6.0	75
G14822	025707-S1	7.2	70
Left Ureter			
G14780	025807	7.2	75
Left Ureter-Open End			
G16817	025806-S1	6.0	75
G14823	025807-S1	7.2	70

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com-A4

