
# Held to a higher standard





## Reach for the single-use solution that offers an enhanced experience from grip to tip.





### IMAGE PROCESSOR

Get enhanced images using the processor's proprietary algorithms, and customize your settings to fine-tune color, saturation, brightness, contrast and more.



### The Ascend Single-Use Flexible Ureteroscope

takes vision versatility to new heights. It's there when and where you need it as your everyday scope, as a backup for when your reusable scope is unavailable, or as your go-to option for those difficult procedures.

### Avoid the pain points associated with reusable scopes:

- Declining functionality (scope fatigue)
- Downtime due to reprocessing
- Cross-contamination risks
- Significant investment and upkeep costs

## **Elevating your endoscopic**

### **expectations.** A complete stone management portfolio is now at your disposal.

Access: Wire Guides

**Dilation:** Sheaths

Lithotripsy: Quanta Lasers

**Extraction:** Nitinol Baskets

Drainage: Catheters and Stents

**Vision:** The Ascend Single-Use Flexible Ureteroscope



### Ascend<sup>™</sup> SINGLE-USE FLEXIBLE URETEROSCOPE

Use the Ascend scope with endoscopic accessories such as laser fibers, stone extractors, biopsy forceps, and the Endoscopic Video Image Processor for endoscopy and endoscopic surgery within the urinary tract and interior of the ureter and kidney.

#### Ascend<sup>™</sup> Single-Use Flexible Ureteroscope

Order Number	Reference Part Number	Working Length mm	Shaft OD Fr	Instrument Channel Fr	Field of View	Deflection Angle	Deflection Configuration
G60337	ASUS-86-300R	670	8.6	3.6	120°	275° up/down	Reverse
G60339	ASUS-86-300S	670	8.6	3.6	120°	275° up/down	Standard

#### **Endoscopic Video Image Processor**

Order Number Reference Part Number		Reference Part Number	LCD screen inch	Dimensions mm	
	G60350	ASUS-IPD-V100	10.1 LCD	264 x 210 x 38	

### **Starter Bundle**

Order Number Reference Part		Reference Part Number	Details
	G60379	ASUS-BUNDLE-10	10 scopes and 1 processor
To order t		То	order the starter bundle, please contact your Cook representative.

The Ascend is not available in all markets. Contact your Cook representative for further information.

As stipulated in our partnership agreement, Cook Medical is only selling Quanta System lasers and fibers in the United States and the following European countries: Austria, France, Germany, Ireland, Switzerland, and the United Kingdom.

#### **Customer Support**

EU Website: cookmedical.eu

EDI: cookmedical.eu/edi Distributors: +353 61239240, ssc.distributors@cookmedical.com Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001702, be.orders@cookmedical.com Denmark: +45 38487607. da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Poland: +48 223060159, pl.orders@cookmedical.com Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

COOK<sup>•</sup> MEDICAL USA Website: cookmedical.com EDI: cookmedical.com/edi.do Americas: Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: customersupport@cookmedical.com Australia: Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283 E-mail: cau.custserv@cookmedical.com

`CC, URO-A4