

Set the standard.



Ascend[®]
SINGLE-USE FLEXIBLE URETEROSCOPE

7.5
FRENCH



Expand your vision portfolio with the **7.5 Fr Ascend Single-Use Flexible Ureteroscope**, and put control in the palm of your hand.

Shaft: 7.5 Fr
outer diameter

285° up-down
deflection




IMAGE PROCESSOR

Get enhanced images using the processor's proprietary algorithms, and customize your settings to fine-tune color, saturation, brightness, contrast, and more.

You're in complete control.

Explore our comprehensive stone management portfolio, and choose from a full complement of devices that accommodate physician preference, from access to drainage.

ILLUMINATION

The distal tip features a unique transparent design, lit with a dual LED, to illuminate the surgical site and the surrounding area.

CUSTOMIZATION

A revamped handle places three action buttons right at your fingertips to tailor in-procedure functionality—such as white-balance adjustment, snapshot, video recording, and more—to your preferred workflow.

DEFLECTION

Static deflection sets your sights on the target.



Vision

Single-use flexible ureteroscopes and cysto-nephroscopes in multiple sizes


Drainage

Catheters and stents




Extraction

Nitinol baskets




Dilation

Sheaths—standard and suction



Access

Wire guides



Lithotripsy

Quanta lasers





Ascend®

SINGLE-USE FLEXIBLE URETEROSCOPE

The Ascend Single-Use Flexible Ureteroscope is developed for use with endoscopic accessories and an Endoscopic Video Image Processor for endoscopy and endoscopic surgery within the urinary tract and the interior of the ureter and kidney, in professional healthcare environments (e.g., hospitals and clinics).

Order Number	Reference Part Number	Working Length mm	Shaft OD Fr	Instrument Channel Fr	Field of View	Deflection Angle	Deflection Configuration
Ascend Single-Use Flexible Ureteroscope							
G61594	ASUS-75-300R	670	7.5	3.6	120°	285° up-down	Reverse
G61595	ASUS-75-300S	670	7.5	3.6	120°	285° up-down	Standard

Please follow any step-by-step instructions carefully, including the risk and technical information provided in the IFU. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.



Order Number	Reference Part Number	Screen	Image Processor Dimensions mm
Endoscopic Video Image Processor HD			
G61562	ASUS-IPD-V200	13.3" LCD	327 x 214 x 54

Please follow any step-by-step instructions carefully, including the risk and technical information provided in the IFU. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Order Number	Reference Part Number	Details
Starter Bundle		
G60379	ASUS-BUNDLE-10	10 scopes and 1 processor

To order the starter bundle, please contact your Cook representative.

Please follow any step-by-step instructions carefully, including the risk and technical information provided in the IFU. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

The Ascend is not available in all markets. Contact your Cook representative for further information.

As stipulated in our partnership agreement, Cook Medical is only selling Quanta System lasers and fibers in the United States and the following European countries: Austria, France, Germany, Ireland, Switzerland, and the United Kingdom.

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 2161 272 3770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do
 Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 Email: customersupport@cookmedical.com
 Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



CC, URO-A4