

Be in control.



Syfonix™
SUCTION URETERAL ACCESS SHEATH

When control and efficiency matter, Syfonix® Suction Ureteral Access Sheath lets you manage the suction at your fingertips during ureteroscopic kidney stone clearance procedures.





Syfonix™

SUCTION URETERAL ACCESS SHEATH

The ureteral access sheath is used to establish a working channel during endoscopic urological procedures, facilitating the passage of endoscopes and other instruments into the urinary tract. The intended population of this device should be adult only. This is a single-use device. A resterilization by users and reuse of the device are not allowed.

Order Number	Reference Part Number	ID Fr	Length cm	Volume mL
G61637	SYUAS-9-40	9	40	-
G61638	SYUAS-9-50	9	50	-
G61639	SYUAS-10-40	10	40	-
G61640	SYUAS-10-50	10	50	-
G61641	SYUAS-11-40	11	40	-
G61642	SYUAS-11-50	11	50	-
G61643	SYUAS-12-40	12	40	-
G61644	SYUAS-12-50	12	50	-

Collection Container

G61646	SYUAS-CC-050	-	-	50
G61647	SYUAS-CC-180	-	-	180

Please follow any step-by-step instructions carefully, including the risk and technical information provided in the IFU. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 2161 272 3770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do

Americas:
 Phone: +1 812.339.2235, 855.207.7214, Fax: 800.330.5490
 Email: customersupport@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



URO-CAP-EQ-8.5X11