

High-definition clarity




Customizable
control buttons
on handle

Ascend[®]

SINGLE-USE FLEXIBLE CYSTO-NEPHROSCOPE HD





Ascend®

SINGLE-USE FLEXIBLE CYSTO-NEPHROSCOPE HD

The Ascend Single-Use Flexible Cysto-Nephroscope HD shall be used in conjunction with the Endoscopic Image Processor HD and other peripheral equipment for observation, diagnosis, photography, and treatment within the bladder, urethra, and kidney.

Order Number	Reference Part Number	Working Length mm	Shaft OD Fr	Instrument Channel Fr	Field of View	Deflection Angle	Deflection Configuration
Ascend Single-Use Flexible Cysto-Nephroscope HD							
G61346	ASUS-HD-15R-CYS	380	15.0	6.6	120°	225° up/down	Reverse
G61347	ASUS-HD-15S-CYS	380	15.0	6.6	120°	225° up/down	Standard

Please follow any step-by-step instructions carefully, including the risk and technical information provided in the IFU. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.



Order Number	Reference Part Number	Screen	Dimensions mm
Endoscopic Video Image Processor HD			
G61562	ASUS-IPD-V200	13.3" LCD	327 x 214 x 54

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 2161 272 3770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do

Americas:
 Phone: +1 812.339.2235, 855.207.7214, Fax: 800.330.5490
 Email: customersupport@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com

URO-CAP-EQ-8.5X11

