




PRESSURE MONITORING CATHETERS



Optimise arterial line placement.



PRESSURE MONITORING CATHETERS



Used for arterial pressure monitoring and blood sampling.

FEATURES & BENEFITS

- Inserted via the minimally invasive Seldinger insertion technique to minimise vessel trauma.
- Lines available for radial or femoral insertion to increase treatment options.
- Femoral access set available with extra-firm nylon catheter for easier insertion through tougher, thicker tissue.

Order Number	Reference Part Number	Catheter Fr/Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Tip Configuration	Needle Baseplate
G02926	C-PMS-250-RA	2.5/21	5	.015	straight	included
G02930	C-PMS-251-RA	2.5/21	5	.015	straight	not included
G01908	C-PMS-300-RA	3.0/20	5	.018	straight	included
G02054	C-PMS-301-RA	3.0/20	5	.018	straight	not included
G02133	C-PMS-301J-RA	3.0/20	5	.018	straight & 2 mm J	not included
G02250	C-PMS-300-FA	3.0/20	8	.018	straight	included
G04021	C-PMS-300J-FA	3.0/20	8	.018	straight & 2 mm J	included
G02598	C-PMS-301-FA	3.0/20	8	.018	straight	not included
G02092	C-PMS-400-RA	4.0/18	6	.021	straight	included
G01909	C-PMS-400-FA	4.0/18	12	.021	straight	included
G01910	C-PMS-400J-FA	4.0/18	12	.021	straight & 2 mm J	included
G02058	C-PMS-401-FA	4.0/18	12	.021	straight	not included
G02089	C-PMS-401J-FA	4.0/18	12	.021	straight & 2 mm J	not included
G02179	C-PMS-500-FA	5.0/16	15	.035	straight	included
G04503	C-PMS-501J-FA	5.0/16	15	.035	straight & 3 mm J	not included
G13332	C-NPMS-501J-15	5.0/16	15	.035	straight & 3 mm J	not included

Sets include arterial catheter, introducer needle, wire guide, dilator (femoral artery access sets only) and syringe.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
 E-mail: cau.custserv@cookmedical.com

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN'S HEALTH

