



Ultraxx™ NEPHROSTOMY BALLOON

The Ultraxx Nephrostomy Balloon features a rated burst pressure of 20 atm and a functional dilating portion 15 cm long. This combination of high pressure and balloon length is designed to establish a nephrostomy tract.



The set includes an Ultraxx balloon catheter, an access sheath, and the Cook Inflation Device.

Features				
Port grip	The grip provides control and aids handling during the inflation process.			
Sheath lock	The lock holds the sheath in place and out of the way during dilation.			
Dual inflation/ deflation ports	The ports provide uniform pressure during inflation. When the balloon is deflated, the ports provide fluid evacuation for the removal process.			
Platinum band	The band provides radiopacity to aid in balloon dilation positioning.			
Tapered tip	The tip accommodates insertion and allows placement over a .038 inch diameter wire guide.			



Ultraxx[™]

The Ultraxx Nephrostomy Balloon is used to dilate the musculofascia, renal capsule, and parenchyma to establish and maintain a percutaneous nephrostomy tract. The set is available with a standard PTFE access sheath or a clear vinyl access sheath for visualisation.

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm		
Set With Clear Vin	yl Sheath				
G46037	UNBS-6-15-CS	6	15		
G46038	UNBS-8-15-CS	8	15		
G30357	UNBS-10-15-CS	10	15		
With PTFE Sheath					
G46036	UNBS-6-15	6	15		
G46024	UNBS-8-15	8	15		
G30365	UNBS-10-15	10	15		

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm		
Balloon Catheter Only With Clear Vinyl Sheath					
G30358	UNB-10-15-CS	10	15		
With PTFE Sheath					
G30364	UNB-10-15	10	15		

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu EDI: cookmedical.eu/edi Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001633, be.orders@cookmedical.com Denmark: +45 38487607. da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Poland: +48 223060159, pl.orders@cookmedical.com Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com EDI: cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: orders@cookmedical.com

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



CC, URO-A4