

Make a diagnosis,
not a repeat biopsy.



BiGopsy[®]
BACKLOADING BIOPSY FORCEPS



cookmedical.eu




BIGopsy®

BACKLOADING BIOPSY FORCEPS

The BIGopsy Backloading Biopsy Forceps are specifically designed to obtain larger sample sizes for the pathological evaluation of tumours within the renal pelvis and ureters. The device may also be used for the retrieval of specimens or foreign objects from the urinary tract. The BIGopsy's 4 mm³ cup size provides a sample larger than any other flexible ureteroscopic cup biopsy forceps on the market that are 3 Fr or smaller. The BIGopsy's larger cup size is enough to increase the likelihood of obtaining a representative tissue sample for the pathologist without costly repeated procedures.

Because of the forceps' superior cup size, BIGopsy has been designed to be backloaded through the working channel of a ureteroscope before being used in the patient. Once the device has been backloaded, secure the removable, replaceable handle to the forceps. The device is then ready for operation.

Note: Deploy this device by using a ureteral access sheath.



Order Number	Reference Part Number	Fr	Length cm
G48240	BLB-024115	2.4	115

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details. Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



©CC, URO-A4