


More options in lithotripsy.



Conquest TTC[®]

LITHOTRIPTOR CABLE

We now offer the Conquest TTC in two flexible cable sizes to give you options for your lithotripsy procedures. And you can use the Conquest with a wide selection of extraction baskets as well as with our Soehendra Lithotriptor Handle. Cook Medical: Giving you the choices you need to help you achieve quality, cost-effective patient care.





Conquest TTC[®]

LITHOTRIPTOR CABLE

Used to mechanically crush stones in the biliary duct when other methods of endoscopic removal have failed.

Conquest is disposable and intended for single use only. The 'through the channel' lithotripter cable accepts baskets compatible with the Soehendra Lithotripter Handle, which is reusable.

- The coil-spring cable provides flexibility and length for use through the channel of flexible endoscopes.
- The flush adapter allows for contrast injection during procedure.
- It is compatible with Soehendra Lithotripter compatible soft wire extraction baskets.

Order Number	Reference Part Number	Fr	Length cm	Minimum Channel Size mm	Comments
G25472	TTCL-1	8.5	170	3.2	Cable
G25518	TTCL-10	10	170	3.7	Cable
G21860	SLH-1	N/A	N/A	N/A	Handle

Soehendra Lithotripter Handle sold separately. Compatible with Memory 7 Fr Soft Wire, The Web[®] Extraction Baskets, and Fusion[®] Lithotripsy Compatible Extraction Basket, sold separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4