

Take complete control of your PEG placement procedures.



Percutaneous Endoscopic Gastrostomy System

From the moment you open the Cook PEG kit, you'll notice the difference. The components are conveniently arranged in the sequential order in which they will be used, making the procedure flow as smoothly and efficiently as the nutrients you are delivering to your patients. The Cook PEG kit—available in push or pull methods—delivers the efficiency you need and the nutrients your patients need.

Percutaneous Endoscopic Gastrostomy System

Used for percutaneous endoscopic placement to provide enteral nutrition to patients requiring nutritional support. It is supplied sterile and is disposable. It is intended for single use only.

- The kit components, including the hemostat and tube clamp, are conveniently packaged for sequential use during the procedure for increased efficiency.
- The clear feeding tube is made of flexible, high-grade medical silicone.
- Reference "X mark" to specify where to cut.
- Compatible, through-the-PEG jejunal tubes are available for all PEG options.

Order Number	Reference Part Number	Feeding Tube Fr	Feeding Tube Inner Diameter mm	Insertion Wire Diameter inch	Insertion Wire Length cm
G22637	FLOW-20-PUSH-S	20	4.47	0.035	260
G22638	FLOW-20-PULL-S	20	4.47	–	150
G22635	PEG-24-PUSH-S	24	5.5	0.035	260
G22636	PEG-24-PULL-S	24	5.5	–	150
G22762	FLOW-20-PUSH-I-S*	20	4.47	0.035	260
G22761	FLOW-20-PULL-I-S*	20	4.47	–	150
G22758	PEG-24-PUSH-I-S*	24	5.5	0.035	260
G22757	PEG-24-PULL-I-S*	24	5.5	–	150

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details. Refer to the current Instructions for Use (IFU) for detailed system use. Please see product risk information in the IFU at cookmedical.eu.

KIT CONTENTS:

1 fenestrated drape	1 silicone feeding tube	1 universal adapter	1 patient care manual	1 ChloraPrep® swabs*
1 150 cm long looped insertion wire (pull only)	5 gauze pads, 4 x 4 inch	1 bolus adapter	1 #11 scalpel	1 povidone ointment*
1 260 cm long wire guide (push only)	1 scissors	2 twist locks	1 trocar needle	1 Xylocaine® 1%, 5.0 mL glass vial*
1 cold snare	1 hemostat	1 cable tie	2 water-soluble lubricant*	1 syringe*
	1 bolster	1 bolus feeding adapter – male (international only)	1 22 gauge x 1½ inch needle*	
	1 tube clamp		1 25 gauge x 5/8 inch needle*	

* International Kits without pharmaceuticals (-I Kits) do not include lubricant, injection needles, syringe, antiseptic swabs, antiseptic ointment, and xylocaine

ChloraPrep is a registered trademark of CareFusion 2000, Inc.

Xylocaine is a registered trademark of Fresenius Kabi USA, LLC.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4