


Access made easy.

Precision delivered
from over 40 years of
proven sphincterotome
technologies.



SureTome™ SW SPHINCTEROTOME WITH DOMETIP®

Building on more than 40 years of proven sphincterotome technology to support endoscopic retrograde cholangiopancreatography (ERCP) procedural success, the SureTome SW Sphincterotome with DomeTip is optimized for short-wire technique, enabling a seamless transition to physician-controlled wire guide cannulation.





SureTome™ SW

SPHINCTEROTOME WITH DOMETIP®

The SureTome SW Sphincterotome with DomeTip is used for the cannulation of the ductal system and for sphincterotomy. If preloaded, also aids in bridging difficult strictures during ERCP.

- Custom 3D forming wire designed for reliable cutting wire orientation
- DomeTip technology designed with anatomy in mind for access
- Protector™ insulated coating may aid in targeting desired tissue to cut
- Rotatable handle with synchronized tip rotation to allow for precise positioning
- Proprietary pre-slit catheter is optimized for use with short-wire technique
- Available as a stand-alone or preloaded with the Acrobat® 2 Calibrated Tip Wire Guide

Order Number	Reference Part Number	Cutting Wire Length mm	Distal Tip Length mm	Distal Tip OD Fr	Catheter OD Fr	Wire Guide Diameter inch	Wire Guide Length cm	Minimal Accessory Channel Diameter mm
Stand-alone Sphincterotome								
G61350	SUR-SW-35	25	7.0	4.4 DomeTip®	7.5	-	-	4.2
Preloaded with the Acrobat® 2 Calibrated Tip Wire Guide*								
G61351	SUR-SW-AWG2-35-260	25	7.0	4.4 DomeTip®	7.5	0.035	260	4.2
G61352	SUR-SW-AWG2-35-450	25	7.0	4.4 DomeTip®	7.5	0.035	450	4.2

**Compatible with wire guide diameter 0.025 and 0.035 inch. Please follow any step-by-step instructions carefully, including the risk and technical information provided in the IFU. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Customer Service

EU website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, emea.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 2161 272 3770, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA website: cookmedical.com
 EDI: cookmedical.com/edi.do

Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 Email: customersupport@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 Email: cau.custserv@cookmedical.com



AI, ESC, IR, PI-8.5X11