Retrieve with ease.



Cook Medical is the only company to offer a 4-loop snare in the market. Indy OTW Vascular Retriever provides a unique over-the-wire capability which allows greater access to vessel coverage, reducing the need for device exchange and saving you valuable time. Other unique features include kink-resistant Flexor® technology and visibility features throughout the entire device.

The Cook difference

Indy C

VASCULAR RETRIEVER

Cook Medical has the industry's broadest line of vascular products. These products and others like them are expanding treatment options and improving patients' lives every day. **Proven products, advancing procedures, clinical partners.**

ACCESS + TARGET	Cook Medical	Medtronic	Gore Medical	Terumo Aortic	Cordis	Endologix	Boston Scientific	AngioDynamics	Merit Medical
Diagnostic Catheters		\bigcirc	\bigcirc	\bigcirc		\bigcirc	\bigcirc		
Wire Guides		\bigcirc	\bigcirc		\bigcirc	\bigcirc			
Molding Balloons				\bigcirc		\bigcirc		\bigcirc	
Sheaths					\bigcirc	\bigcirc		\bigcirc	
Selective Catheters		\bigcirc	\bigcirc			\bigcirc	\bigcirc		
TREAT	Cook Medical	Medtronic	Gore Medical	Terumo Aortic	Cordis	Endologix	Boston Scientific	AngioDynamics	Merit Medical
Abdominal							\bigcirc	\bigcirc	\bigcirc
Thoracic					\bigcirc	\bigcirc	\bigcirc	\bigcirc	0
Visceral		\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc



Competitive analysis performed by Cook Medical (February 2024). This chart is for illustration purposes.

PRODUCT RELEASE

Vendor name	Cook Medical		Date		
Contract					
Contract catego	ory		Contract number		
Product name	Indy OTW Vascular Retrie	ever	Product number		
Contract price			Unit of sale 1		
Service line					
Product description					
The Indy OTW Vascular Retriever consists of a 4-loop wire snare attached to an inner catheter (which is attached proximally to a metal cannula) and an outer Flexor sheath. The device is compatible with a 0.035 inch wire guide for over-the-wire introduction. The wire snare is composed of nitinol for shape memory, and platinum for visibility. The Flexor sheath has a radiopaque band incorporated within the sheath material to identify the location of the sheath's distal tip.					

Product is: Reusa	ible Disposable	Shell Life:	Latex Free: Yes INO
ITEMS TO BE RE	PLACED		
Ve	ndor	Product Name	Product Description
Desite and at a second	lean a ta alata na mandura a c	and Iliac branch devices fenestrated complex an	rtic or general procedures where a

During what procedures is this product used? Illac branch devices, fenestrated, complex aortic or general procedures w snare is needed to retrieve a foreign body

What advanced features does this product offer over current methods?

Over-the-wire introduction, Flexor technology, 4-loop design

What supply chain advantages does this product offer over current methods?

Cook is the only company to offer a full array of devices needed to support endovascular procedures. Streamline your purchasing with one reliable supplier.

COOK MEDICAL INCORPORATED P.O. Box 4195, Bloomington, IN 47402-4195 USA Phone: 812.339.2235, Toll Free: 800.457.4500, Toll-Free Fax: 800.554.8335

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



Scan the code for more product information:

