# Perform two procedures with one device.




Illustration by Lisa Clark


# Goldstein SONOBIOPSYTT CATHETER





### Improve the visualization of the uterine cavity and obtain an endometrial biopsy.

The Goldstein SonoBiopsy Catheter allows physicians to perform both sonohysterography and endometrial biopsy in one procedure. This specially designed catheter has a round, closed tip that may make it easier to access the uterus, and has a wide sideport that facilitates the transition between the instillation of saline solution and the aspiration of endometrial tissue.




Illustrations by Lisa Clark

### Sonohysterography

Wide sideport facilitates saline instillation.

# Goldstein

The Goldstein SonoBiopsy Catheter is conveniently marked in 1 cm gradations up to 10 cm to indicate proper positioning. The catheter also features an adjustable acorn-shaped positioner that fits snugly into the cervical os, minimizing leakage.





## Goldstein


#### SONOBIOPSY™ CATHETER

Used to access the uterine cavity for sonohysterography and to obtain an endometrial biopsy, if indicated, using the same device.

The adjustable acorn-shaped positioner maintains catheter placement while reducing leakage of instilled saline at the external cervical os. Ink marks in 1 cm gradations up to 10 cm indicate catheter positioning.

The catheter's 7.2 French size allows it to pass into the uterus for tissue collection. The catheter's round, closed tip eases cervical passage. The large oval sideport facilitates instillation of saline solution and aspiration of endometrial tissue.

Visit www.cookmedical.com/rh to view a step-by-step procedural video.



Order	Reference	Catheter	Length	Quantity
Number	Part Number	Fr	cm	
G46190	J-GSBX-072026	7.2	26	10

 $Some\ products\ or\ part\ numbers\ may\ not\ be\ available\ in\ all\ markets.\ Contact\ your\ local\ Cook\ representative\ or\ Customer\ Service\ for\ details.$ 

### **Customer Service**

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

**Americas:** EDI - www.cookmedical.com/edi.do Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: orders@cookmedical.com

#### Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283 E-mail: cau.custserv@cookmedical.com

