




Zilver 635[®]

BILIARY SELF-EXPANDING STENT




Trust a dedicated leader in drainage products,
featuring the Zilver 635[®] biliary stent.

The Zilver 635 Biliary Self-Expanding Stent has been designed to provide high radial force. The 40 cm delivery system is perfectly suited for percutaneous approaches to the liver. Cook Medical's complete offering of accessory products provides support of your diverse procedural needs.

A precise and durable stent, designed for better results for your patients.

The 6.0 Fr Zilver balances the most important stent traits: foreshortening resistance, radial force, flexibility and durability. You can rely on this stent for precise delivery and excellent patency rates,¹ providing you with confidence in the stent's performance and allowing you to focus on providing better results for your patient.




PRECISION

No clinically significant foreshortening during deployment. Gold markers enhance visibility during fluoroscopy, and electropolishing results in a smooth and durable stent surface.



STRENGTH AND FLEXIBILITY

Z-stent cell design ensures optimal wall apposition while maintaining a strong, equally distributed radial force throughout the stent, and provides excellent kink-resistant attributes within tortuous anatomy without lumen reduction.



RELIABILITY

Self-expanding nitinol construction provides the flexibility and high radial force needed in compressed and nonlinear biliary ducts.


1. Loew BJ, Howell DA, Sanders MK, et al. Comparative performance of uncoated, self-expanding metal biliary stents of different designs in 2 diameters: final results of an international multicenter, randomized, controlled trial. *Gastrointest Endosc.* 2009;70(3):445-453.

A flexible and efficient delivery system, designed specifically for you.

The Zilver 635 Biliary Self-Expanding Stent and its delivery system are designed for ease of use to make your biliary drainage procedures efficient and effective.

- Flexor® technology gives the delivery system the flexibility, kink resistance and pushability needed in tortuous biliary anatomy.
- Flexible 40 cm, 6.0 Fr delivery system is specifically designed for percutaneous biliary procedures.
- A convenient pullback handle and radiopaque marker at the distal tip allow precise and controlled deployment of the stent for smooth delivery into the target anatomy.


◀ Available with a convenient 40 cm Flexor delivery system.



Flexor sheath


A full complement of accessory products.

For a complete biliary procedure, you need devices that help you access, target and treat. With our full range of biliary drainage products, Cook Medical is your procedural partner from start to finish.



ACCESS

Enable minimally invasive access to the biliary tree with needles and introducers designed for precision and smooth access.



TARGET

Achieve accurate targeting and crossing of obstructions with versatile wire guides and catheters.



TREAT

Complete your biliary procedures with a drainage catheter featuring supersmooth Ultrathane® material for patient comfort, or a nitinol stent designed for precision, radial strength and visibility.





Zilver 635[®]

BILIARY SELF-EXPANDING STENT

Order Number	Reference Part Number	Stent Diameter mm	Stent Length cm
G37269	ZIB6-40-8-4.0	8	4
G37270	ZIB6-40-8-6.0	8	6
G37271	ZIB6-40-8-8.0	8	8
G37272	ZIB6-40-9-4.0	9	4
G37273	ZIB6-40-9-6.0	9	6
G37274	ZIB6-40-9-8.0	9	8
G37275	ZIB6-40-10-4.0	10	4
G37276	ZIB6-40-10-6.0	10	6
G37277	ZIB6-40-10-8.0	10	8
G39719	ZIB6-40-12-4.0	12	4
G39720	ZIB6-40-12-6.0	12	6
G39721	ZIB6-40-12-8.0	12	8

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH



Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com
www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

