
ERCP is all about smooth, easy access.



Fusion® OMNI™

SPHINCTEROTOME WITH DOMETIP®

DUAL PLATFORM ERCP | ACCESS | STONE MANAGEMENT | STRICTURE MANAGEMENT

COOK® MEDICAL

Image courtesy of Dr. David Diehl,
Geisinger Health System, Danville, PA.



Fusion® OMNI™

SPHINCTEROTOME WITH DOMETIP®

Used for cannulation of the ductal system and for sphincterotomy.

- DomeTip designed with anatomy in mind for access
- Custom 3D forming wire designed for reliable cutting wire orientation
- Breakthrough channel technology provides options for the physician to control the wire guide or use the device with an over-the-wire long wire technique
- Preloaded options are available for procedural convenience with a 0.021 or 0.035 inch diameter wire guide

Order Number	Reference Part Number	Catheter Fr	Tip Shape	Wire Guide Diameter inch	Wire Guide Length cm	Minimum Accessory Channel Diameter mm
G31903	FS-OMNI*	7.5	DomeTip	0.035	–	3.2
G48675	FS-OMNI-21**	5.4	DomeTip	0.021	–	2.8
Preloaded with Acrobat®2 Calibrated Tip Wire Guide						
G57301	FS-OMNI-AWG2-35-205	7.5	DomeTip	0.035	205	3.2
G57302	FS-OMNI-AWG2-35-260	7.5	DomeTip	0.035	260	3.2
G57303	FS-OMNI-AWG2-35-450	7.5	DomeTip	0.035	450	3.2
Preloaded with Tracer Metro® Direct™ Wire Guide						
G56103	FS-OMNI-21-260	5.4	DomeTip	0.021	260	2.8
G48676	FS-OMNI-21-480	5.4	DomeTip	0.021	480	2.8
G55814	FS-OMNI-35-260	7.5	DomeTip	0.035	260	3.2
G55815	FS-OMNI-35-480	7.5	DomeTip	0.035	480	3.2

Cutting wire length is 25 mm.
**Recommended:Acrobat® 2 Calibrated Tip Wire Guide, sold separately. **Recommended: Tracer Metro® Direct™ Wire Guide, sold separately.*

Please follow any step-by-step instructions carefully, including the risk, technical information, Device Name and intended use provided in the IFU at cookmedical.eu. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
EDI: cookmedical.com/edi.do

Americas:
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: customersupport@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

COOK® MEDICAL

AI-ESC-IR-OHNS-PI-RH-SUR-A4

© COOK 07/2026 ESC-WF1254700-EN