
An alternative for ductal access.




EchoTip® Ultra

HD ULTRASOUND ACCESS NEEDLE

When you have a pancreatobiliary access case that is too challenging for ERCP, you need an alternative. That's why Cook Medical created the HD Ultrasound Access Needle—the first therapeutic EUS needle indicated for gaining access. And, because it has our highly echogenic dimple pattern and ergonomic handle, it can help assure that the needle is within the target and that the procedure is stable and secure.

COOK® MEDICAL



EchoTip® Ultra

HD ULTRASOUND ACCESS NEEDLE

This device is used to access or sample submucosal and extramural lesions of the gastrointestinal tract, access of the following: the intra- or extra-hepatic bile ducts, pancreatic ducts, cystic duct, gallbladder or for delivery of injectable materials into tissues through the accessory channel of an ultrasound endoscope.

Supplied sterile and is disposable – intended for single use only.

- Dimple pattern promotes ultrasound visibility
- Naturally contoured handle delivers precise needle control and stability
- Atraumatic needle tip designed to allow passage of a .035" wire guide to maintain access
- Beveled stylet used to puncture and gain access
- Stylet cap designed to maintain stylet position on puncture

Order Number	Reference Part Number	Needle Gage	Needle Extension cm	Sheath Fr	Stylet Tip	Minimum Accessory Channel mm
G52012	ECHO-HD-19-A	19	0-8	5.2- 4.3	Beveled	2.8

5 cc and 10 cc vacuum locking syringe included.

Please follow any step-by-step instructions carefully, including the risk, technical information, Device Name and intended use provided in the IFU at cookmedical.eu. Note that some products or part numbers may not be available in all markets; consult your Cook representative for details.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, is.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland – French: +41 448009609, fr.orders@cookmedical.com
Switzerland – Italian: +41 448009609, it.orders@cookmedical.com
Switzerland – German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
EDI: cookmedical.com/edi.do
Americas:
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com
Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

AI-ESC-IR-OHNS-PI-RH-SUR-A4

© COOK 04/2026